

Product Datasheet

Cytomegalovirus (CMV-ICP36) Antibody CA006

Product Name	Cytomegalovirus (CMV-ICP36) Antibody
Catalog Number	CA006
Species/Host	Mouse
Tested Applications	ELISA, WB, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	ICP36