

Product Datasheet

Cytomegalovirus (CMV) Antigen, Native CV046

Product Name	Cytomegalovirus (CMV) Antigen, Native
Catalog Number	CV046
Tested Applications	WB, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤-70°C
Note	Confirmed negative for HIV and Hepatitis B Virus.
Purity	Density gradient purified
Source	NHDF Cells/Strain AD169
Buffer	10 mM Tris, 100 mM NaCl, 1 mM EDTA, 0.1% Triton X-100, pH 8