

Product Datasheet

Measles Virus Antigen, Native EK501

Product Name	Measles Virus Antigen, Native
Catalog Number	EK501
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Grade 2 Antigen
Storage	≤-70°C. Avoid repeated freezing and thawing.
Note	Sonicate immediately prior to usage. Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Purity	Partially purified
Source	Vero Cells/Strain Edmonston
Buffer	MEM