

Product Datasheet

Murine IP-10 ELISA Kit 660.120

Product Name	Murine IP-10 ELISA Kit
Catalog Number	660.120
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	7.8-500pg/ml
Sensitivity	6.5pg/ml
Tested Applications	ELISA
Alternative Names	CXCL10
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IP-10
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable