

Product Datasheet

Thyroid Stimulating Hormone β (TSH β) Antibody FA305

Product Name	Thyroid Stimulating Hormone β (TSH β) Antibody
Catalog Number	FA305
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Human TSH β -subunit
Alternative Names	Thyrotropin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 %
Buffer	0.05M Phosphate Buffered Saline, pH 7.2 \pm 0.2
Cross Reactivity	TSH 100%, β hCG<1%, ihCG<1%, LH<1%, FSH<1%
Specificity	TSH β
Matched Pairs	FA316 (detection)
Preservatives	0.1% Sodium Azide