

## Product Datasheet

### Luteinizing Hormone $\beta$ (LH- $\beta$ ) Antibody FA309

<b>Product Name</b>	Luteinizing Hormone $\beta$ (LH- $\beta$ ) Antibody
<b>Catalog Number</b>	FA309
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Human Luteinizing hormone $\beta$ -subunit
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	0.05M Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	LH 100%, $\beta$ hCG > 10%, ihCG >1%, FSH < 1%, TSH < 1%
<b>Specificity</b>	Human Luteinizing hormone $\beta$ -subunit
<b>Matched Pairs</b>	FA313 (detection)
<b>Preservatives</b>	0.1% Sodium Azide