

Product Datasheet

Luteinizing Hormone β (LH- β) Antibody FA309

Product Name	Luteinizing Hormone β (LH- β) Antibody
Catalog Number	FA309
Description	Purified mouse monoclonal antibody
Species/Host	Mouse
Tested Applications	ELISA
Alternative Names	LH, Lutropin, Interstitial Cell-Stimulating Hormone
Brand	EastCoast Bio
Storage	+2-8 °C
Isotype	Mouse IgG
Clonality	Monoclonal
Buffer	0.05 M phosphate-buffered saline, pH 7.2 \pm 0.2, 0.1% sodium azide as a preservative.
Cross Reactivity	LH 100%, β hCG > 10%, ihCG >1%, FSH < 1%, TSH < 1%
Specificity	Beta subunit of human luteinizing hormone (lutropin, LH β)
Matched Pairs	HM222 (detection)