

Product Datasheet

Prolactin (PRL) Antibody FA310

Product Name	Prolactin (PRL) Antibody
Catalog Number	FA310
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Human Prolactin
Alternative Names	Prolactin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG
Clonality	Monoclonal
Buffer	0.05M Phosphate Buffered Saline, pH 7.2 ± 0.2
Specificity	Human Prolactin
Matched Pairs	HM120 (capture)
Preservatives	0.1% Sodium Azide