

## Product Datasheet

### Prolactin (PRL) Antibody FA312

<b>Product Name</b>	Prolactin (PRL) Antibody
<b>Catalog Number</b>	FA312
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Human Prolactin
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	0.05M Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	Human Prolactin
<b>Matched Pairs</b>	HM120 (capture)
<b>Preservatives</b>	0.1% Sodium Azide