

Product Datasheet

Triiodothyronine (T3) Antibody FA315

Product Name	Triiodothyronine (T3) Antibody
Catalog Number	FA315
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	T3-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C. Do not freeze.
Isotype	IgG
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.05M Phosphate Buffered Saline, pH 7.2 ± 0.2
Specificity	T3
Preservatives	0.1% Sodium Azide