

## Product Datasheet

### Human Immunoglobulin E (IgE) Antibody FA317

<b>Product Name</b>	Human Immunoglobulin E (IgE) Antibody
<b>Catalog Number</b>	FA317
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Human IgE
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	+2-8 °C. Do not freeze.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 %
<b>Buffer</b>	0.05M Phosphate Buffered Saline, pH 7.2 ± 0.2
<b>Specificity</b>	Human IgE
<b>Matched Pairs</b>	FA318 (capture)
<b>Preservatives</b>	0.1% Sodium Azide