

Product Datasheet

Human Immunoglobulin E (IgE) Antibody FA317

Product Name Human Immunoglobulin E (IgE) Antibody

Catalog Number FA317

Species/Host Mouse

Tested Applications ELISA

Immunogen Human IgE

Brand EastCoast Bio

Form/Appearance Purified Monoclonal Antibody

Storage 2-8°C. Do not freeze.

Isotype IgG2b

Clonality Monoclonal

Purity >95%

Buffer 0.05 M Phosphate Buffered Saline, pH 7.4

Specificity Human IgE

Matched Pairs FA318 (capture)

Preservatives 0.1% Sodium Azide