

Product Datasheet

Human Immunoglobulin E (IgE) Antibody FA318

Product Name Human Immunoglobulin E (IgE) Antibody

Catalog Number FA318

Species/Host Mouse

Tested Applications ELISA

Immunogen Human IgE

Brand EastCoast Bio

Form/Appearance Purified Monoclonal Antibody

Storage 2-8°C

Isotype IgG2b

Clonality Monoclonal

Purity >95%

Buffer 0.05 M Phosphate Buffered Saline, pH 7.4

Specificity Human IgE

Matched Pairs FA317 (detection)

Preservatives 0.1% Sodium Azide