

Product Datasheet

Toxoplasma gondii Antigen, Native FC511

| | |
|---------------------|--|
| Product Name | Toxoplasma gondii Antigen, Native |
| Catalog Number | FC511 |
| Tested Applications | ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid, sonicated and centrifuged tachyzoites |
| Storage | ≤-70°C |
| Note | Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious. |
| Source | RH Strain in PBS |
| Buffer | 10 mM PBS, pH 7.4 with 1 mM AEBSF and 10 KIU/ml Aprotinin |