

Product Datasheet

Toxoplasma gondii Antigen, Native FC511

Product Name	Toxoplasma gondii Antigen, Native
Catalog Number	FC511
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid, sonicated and centrifuged tachyzoites
Storage	≤-70°C
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	RH Strain in PBS
Buffer	10 mM PBS, pH 7.4 with 1 mM AEBSF and 10 KIU/ml Aprotinin