

# Product Datasheet

## Rat IL-6 ELISA Kit 670.010

Product Name	Rat IL-6 ELISA Kit
Catalog Number	670.010
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.3-2000 pg/ml
Sensitivity	12 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat IL-6
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable.