

## Product Datasheet

### Adenovirus Antibody FG300

<b>Product Name</b>	Adenovirus Antibody
<b>Catalog Number</b>	FG300
<b>Description</b>	Protein G Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	WB, ELISA, LF, IFA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein G Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0-7.5
<b>Specificity</b>	Hexon Protein
<b>Matched Pairs</b>	HM1178 (detection)