

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody FK100

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	FK100
Description	Immunoaffinity Purified Polyclonal Antibody,
Species/Host	Goat
Tested Applications	ELISA, LF
Immunogen	HBsAg, ad and ay
Brand	EastCoast Bio
Form/Appearance	Immunoaffinity Purified Polyclonal Antibody,
Storage	+2-8 °C
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	HBsAg, ad and ay
Matched Pairs	HM301 (detection)
Preservatives	0.02% Sodium Azide