

Product Datasheet

Influenza A Virus Antibody FT351

Product Name	Influenza A Virus Antibody
Catalog Number	FT351
Species/Host	Mouse
Tested Applications	ELISA, WB, IHC
Immunogen	Influenza A Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Cross Reactivity	No reactivity with Influenza B, Adenovirus, and RSV
Specificity	Nucleoprotein from all strains
Preservatives	0.02% Sodium Azide