

Product Datasheet

Influenza A Virus Antibody FT351

| | |
|---------------------|---|
| Product Name | Influenza A Virus Antibody |
| Catalog Number | FT351 |
| Species/Host | Mouse |
| Tested Applications | ELISA, WB, IHC |
| Immunogen | Influenza A Virus |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | >95% by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.0 |
| Cross Reactivity | No reactivity with Influenza B, Adenovirus, and RSV |
| Specificity | Nucleoprotein from all strains |
| Preservatives | 0.02% Sodium Azide |