

## Product Datasheet

### Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody FT372

<b>Product Name</b>	Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody
<b>Catalog Number</b>	FT372
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, IHC
<b>Immunogen</b>	Human milk fat globules
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2b, kappa
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffer, pH 7.0
<b>Specificity</b>	MUC1
<b>Preservatives</b>	0.02% Sodium Azide