

# Product Datasheet

## Rat TGF beta 1 ELISA Kit 670.020

Product Name	Rat TGF beta 1 ELISA Kit
Catalog Number	670.020
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	62.5 à 2000 pg/ml
Sensitivity	48 pg/ml
Tested Applications	ELISA
Alternative Names	TGFb1, TGF-beta1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat TGF-b1
Sample Types	Serum, Plasma, Cell culture supernatant