

# Product Datasheet

## Rat IL-1 beta ELISA Kit 670.040

Product Name	Rat IL-1 beta ELISA Kit
Catalog Number	670.040
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.2-2000 pg/ml
Sensitivity	4.4 pg/ml
Tested Applications	ELISA
Alternative Names	IL-1b, IL-1beta
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes both natural and recombinant rat IL-1b
Sample Types	Serum, Plasma, Cell culture supernatant