

# Product Datasheet

## Rat IL-10 ELISA Kit 670.070

Product Name	Rat IL-10 ELISA Kit
Catalog Number	670.070
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000pg/ml
Sensitivity	1.5pg /ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant rat IL-10
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable