

Product Datasheet

Amphetamine Antibody GN300

Product Name	Amphetamine Antibody
Catalog Number	GN300
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 \pm 0.2
Specificity	Amphetamine
Preservatives	0.05% Sodium Azide