

Product Datasheet

Rat IL-17A ELISA Kit 670.080

Product Name	Rat IL-17A ELISA Kit
Catalog Number	670.080
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	1.6-100pg/ml
Sensitivity	1.0pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant rat IL-17A
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable