

## Product Datasheet

### Rat IL-17A ELISA Kit 670.080

<b>Product Name</b>	Rat IL-17A ELISA Kit
<b>Catalog Number</b>	670.080
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	1.6-100pg/ml
<b>Sensitivity</b>	1.0pg/ml
<b>Tested Applications</b>	ELISA
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Specificity</b>	Recognizes both natural and recombinant rat IL-17A
<b>Sample Types</b>	Cell culture supernatant and serum tested, other body fluids suitable