

Product Datasheet

Herpes Simplex Virus 1 (HSV 1-gE) Antibody H1A054

Product Name	Herpes Simplex Virus 1 (HSV 1-gE) Antibody
Catalog Number	H1A054
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gE