

Product Datasheet

Rat IL-1 alpha ELISA Kit 670.090.192

Product Name	Rat IL-1 alpha ELISA Kit
Catalog Number	670.090.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	15.6-1000pg/ml
Sensitivity	4pg/ml
Tested Applications	ELISA
Alternative Names	IL-1a
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant rat IL-1 α
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable