

Product Datasheet

Rat IL-2 ELISA Kit 670.100

Product Name	Rat IL-2 ELISA Kit
Catalog Number	670.100
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	62.5-2000pg/ml
Sensitivity	11pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant rat IL-2
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable