

Product Datasheet

Rat VEGF-C ELISA Kit 670.110

Product Name	Rat VEGF-C ELISA Kit
Catalog Number	670.110
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	47-3000pg/ml
Sensitivity	27pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant rat VEGF-C
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable