

Product Datasheet

Rat VEGF-C ELISA Kit 670.110

| | |
|---------------------|--|
| Product Name | Rat VEGF-C ELISA Kit |
| Catalog Number | 670.110 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 47-3000pg/ml |
| Sensitivity | 27pg/ml |
| Tested Applications | ELISA |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Specificity | Recognizes both natural and recombinant rat VEGF-C |
| Sample Types | Cell culture supernatant and serum tested, other body fluids suitable |