

Product Datasheet

Herpes Simplex Virus 1 & 2 (HSV 1 & 2-gD) Antibody HA025

Product Name	Herpes Simplex Virus 1 & 2 (HSV 1 & 2-gD) Antibody
Catalog Number	HA025
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gD of HSV-1 & 2