

## Product Datasheet

### Ketamine Antibody HM006

<b>Product Name</b>	Ketamine Antibody
<b>Catalog Number</b>	HM006
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF
<b>Immunogen</b>	Ketamine-BSA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	DEAE Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Potassium Phosphate, 0.85% NaCl, pH 7.2
<b>Cross Reactivity</b>	Ketamine 100% NorKetamine: 60%
<b>Specificity</b>	Ketamine
<b>Matched Pairs</b>	LA010 (capture)
<b>Preservatives</b>	0.05% Sodium Azide