

## Product Datasheet

### Methamphetamine Antibody HM008

**Product Name** Methamphetamine Antibody

**Catalog Number** HM008

**Species/Host** Mouse

**Tested Applications** LF

**Brand** EastCoast Bio

**Form/Appearance** Protein A Purified Monoclonal Antibody

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Isotype** IgG

**Clonality** Monoclonal

**Buffer** 10 mM Phosphate Buffered Saline, pH 7.4

**Specificity** Methamphetamine

**Matched Pairs** LA013 or LA014 (capture)

**Preservatives** 0.1% Sodium Azide