

Product Datasheet

Giardia lamblia (G. lamblia cysts) Antibody HM014

Product Name	Giardia lamblia (G. lamblia cysts) Antibody
Catalog Number	HM014
Species/Host	Mouse
Tested Applications	LF, ELISA, WB, IFA
Immunogen	Native G. lamblia cysts
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	65 kDa Giardia antigen
Matched Pairs	HM013, HM015 (detection)
Preservatives	0.1% Sodium Azide