

# Product Datasheet

## Human IL-1 alpha ELISA Kit 850.005

Product Name	Human IL-1 alpha ELISA Kit
Catalog Number	850.005
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	10 pg/ml
Tested Applications	ELISA
Alternative Names	IL-1a
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-1a
Sample Types	Serum, Plasma, Cell culture supernatant