

Product Datasheet

Clostridium difficile toxin B Antibody HM020

Product Name	Clostridium difficile toxin B Antibody
Catalog Number	HM020
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Minimal cross reactivity with toxin A.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	C. difficile toxin B
Matched Pairs	HM021 (detection)
Preservatives	0.1% Sodium Azide