

Product Datasheet

Norovirus genogroup I Antibody HM022

Product Name	Norovirus genogroup I Antibody
Catalog Number	HM022
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Does not cross-react with Norovirus Group II.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Group I capsid protein
Matched Pairs	HM023 (capture)
Preservatives	0.1% Sodium Azide