

Product Datasheet

Norovirus genogroup II Antibody HM025

Product Name	Norovirus genogroup II Antibody
Catalog Number	HM025
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Possible cross-reaction with Norovirus Group I.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Group II capsid protein
Matched Pairs	HM024 (detection)
Preservatives	0.1% Sodium Azide