

Product Datasheet

Dengue Virus (DENV-pan NS1) Antibody HM026

Product Name	Dengue Virus (DENV-pan NS1) Antibody
Catalog Number	HM026
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Immunogen	Recombinant Dengue NS1
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	NS1 protein of all 4 serotypes
Matched Pairs	HM027, HM164 or HM538 (detection)
Preservatives	0.1% Sodium Azide