

Product Datasheet

Dengue Virus (DENV-pan NS1) Antibody HM027

Product Name	Dengue Virus (DENV-pan NS1) Antibody
Catalog Number	HM027
Species/Host	Mouse
Tested Applications	ELISA, IFA
Immunogen	Recombinant Dengue NS1
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	NS1 protein of all 4 serotypes
Matched Pairs	HM026 (capture)
Preservatives	0.1% Sodium Azide