

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody HM029

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody
Catalog Number	HM029
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	HIV-1 p24
Matched Pairs	HM030 (detection)
Preservatives	0.1% Sodium Azide