

Product Datasheet

Chlamydia trachomatis (LPS) Antibody HM031

Product Name	Chlamydia trachomatis (LPS) Antibody
Catalog Number	HM031
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	LPS of Chlamydia Genus (C. trachomatis, C. pneumoniae and C. psittaci)
Matched Pairs	HM031 (self-pairs), HM032 or HM215 (capture)
Preservatives	0.1% Sodium Azide