

Product Datasheet

Campylobacter jejuni Antibody HM035

Product Name	Campylobacter jejuni Antibody
Catalog Number	HM035
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Non-reactive with : H. pylori, E. coli, Salmonellae, Shigellae, P. aeruginosa, Yersinia and Citrobacter
Specificity	Whole Cell
Matched Pairs	HM035 (self-pairs)
Preservatives	0.1% Sodium Azide