

Product Datasheet

Salmonella LPS Core Antibody HM038

Product Name	Salmonella LPS Core Antibody
Catalog Number	HM038
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	No cross-reactivity with E. coli, Klebsiella, Citrobacter, Pseudomonas, Yersinia, Shigella, Proteus or Legionella.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with related bacteria
Specificity	LPS core
Matched Pairs	HM039 or HM040 (detection)
Preservatives	0.1% Sodium Azide