

Product Datasheet

Metapneumovirus Antibody HM043

Product Name	Metapneumovirus Antibody
Catalog Number	HM043
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No cross-reactivity with RSV, Influenza A, Influenza B, Adenovirus, Parainfluenza 1,2 &3
Specificity	Matrix Protein (genotypes A and B)
Matched Pairs	HM044 (detection)
Preservatives	0.1% Sodium Azide