

Product Datasheet

Mumps Virus Antibody HM050

Product Name	Mumps Virus Antibody
Catalog Number	HM050
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No cross-reactivity with Measles Virus, RSV, Parainfluenza 1, 2, 3
Specificity	Mumps Virus nucleoprotein
Preservatives	0.1% Sodium Azide