

Product Datasheet

Mumps Virus Antibody HM051

Product Name Mumps Virus Antibody

Catalog Number HM051

Species/Host Mouse

Tested Applications LF, ELISA, IFA

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG2b

Clonality Monoclonal

Purity >90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Cross Reactivity No cross-reactivity with Measles Virus, RSV, Parainfluenza 1, 2, 3

Specificity Mumps Virus nucleoprotein

Preservatives 0.1% Sodium Azide