

Product Datasheet

Adenovirus Antibody HM056

Product Name	Adenovirus Antibody
Catalog Number	HM056
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified human adenovirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 95 %
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	Adenovirus
Matched Pairs	HM057 (detection)
Preservatives	0.02% Sodium Azide